


Shipping samples using used tip boxes

When shipping samples to us, we recommend using **0.2 ml 8-tube strips**, mainly because it minimizes plastic waste and ethanol use. In addition 0.2 ml tubes tend to seal better than larger tubes, thus minimizing the chance of an ethanol leak.

A **used pipette tip box** provides an easy, and free, way to keep your tube strips organized (see below). We recommend **taping** the tube strips down so that they do not fall out during shipment.



Don't forget to **label your strips**. For example you can number your tube strips 1 through 12 as shown above.

When your samples have been organized as described above, all you need to do is package your tip box and ship it to us using your preferred shipping method (USPS, FedEx, UPS, etc).

For any questions contact us at 1-877-743-6689 or by email at cs@mousegenotype.com.